

DIE CUTTER FOR RING SAMPLES OF CONDOMS AND GLOVES

The Die Cutter for ring samples is designed for rapid preparation of faultless ring samples of condoms and gloves.

The Die Cutter consists of a holder with integrated air pump, stamping sheet and folding blade holder for standard cutters.



Preparation of samples

- Place or transfer a smooth part of paper in the cutting area on the stamping plate,
- mount the condom on the condom holder,
- inflate the condom with the plunger,
- place the cutting die in its cradle,
- move the die cutter in the cutting press.

The Die Cutter enables a quickly preparation of faultless ring samples without previous removal of lubricant or powder.

The inflation with air will prevent the sample from forming wrinkles at cutting.

The Die Cutter will highly reduce the time-consuming adjustment of the specimen.

Specifications

Length:	350 mm
Width:	250 mm
Height:	60 mm
Weight:	1,2 kg

Art. No.: 801

VALENDOR reserves the right to modify this specification in part or as a whole